

## HS-010 Head Storage Box for CCU-010 accessories



The HS-010 vacuum storage box allows you to keep a second coating head for your CCU-010 Compact Coating Unit, an additional planetary or rotary stage and the full range of sputtering targets and carbon thread always clean and under vacuum. An integrated pump line attachable to the basic unit pumps down the storage box to vacuum using the available pumping system of the CCU-010. Alternatively one can easily connect an external vacuum pump instead. A comfortable manual lock and vent valve allows the independent pump and vent control of the storage box. The additional standard vacuum chamber helps to completely separate the carbon coating from metal deposition which effectively avoids cross contamination.





**Technical data**

Depth x With x Height	38 x 12 x 19 cm
Weight	4.5 kg
Number of targets to store	6
Stages to store	1 x PS-010 or 1 x RS-010 or 2 x fix platform
Pump line connection	outer Ø 10 mm / DN 100 ISO-K
Compatibility	SP-010, SP-011, CT-010 CT-010 with option glow discharge

